

Touchless –
a real affinity with hands



 **DÜRR
DENTAL**
THE BEST, BY DESIGN

Touchless universal dispenser

The failsafe Touchless 1000 (for 800 ml bottles) and Touchless 400 (for 500 ml bottles) sensor dispensers have several advantages – functionality, reliability and durability. This is important as hands-free handling prevents cross-contamination. The dispensers provide everything required to ensure smooth procedures in a dental practice.

Touchless T 1000/400 dispensers have been developed from the tried-and-tested universal dispensers with a lever for manual handling. With its failsafe sensor technology, the touchless version retains the advantages of the manual dispenser but with the following added benefits:

- Optimum hand hygiene thanks to contactless operation
- Battery operation with battery power indicator by means of an audible battery replacement signal
- User-friendly dosing technology





Autoclavable pump

Easy replacement from the front.

Flexible sealing cap

Ensures that the bottle sits firmly in the dispenser.

High standards of hygiene

The case is made from anodised aluminium - its bacteriostatic effect hinders the growth of bacteria.

Technical data

Battery voltage 7.5 V DC

Battery type (5x) mono cell ANSI D (LR 20/1.5 V)

Battery service life up to 2 years

Dosage (adjustable):

0.75 ml/1 ml/1.25 ml/1.5 ml

Range of detection of fail-safe sensor:

4 cm from top and bottom

Weight without batteries 1250 g

Dimensions W x H x D: 9.2 x 33.0 x 21.0 cm

Accessories



Drip tray

With stainless steel bracket and removable plastic dish.



Table stand


Made of sturdy plastic for the flexible use of the Touchless universal dispenser.



HD products from Dür System Hygiene

- HD 410 Hand disinfection (500 ml or 800 ml)
- HD 435 Washing lotion (500 ml or 1 l)
- HD 440 Care lotion (500 ml)

Viruses destroyed? Yes!
Skin damaged? No!
Hand disinfection
must be gentle.



To the usage video!

DÜRR DENTAL SE
Höpfigheimer Strasse 17
74321 Bietigheim-Bissingen
Germany
www.duerrdental.com
info@duerr.de

