


## Benefits of the Dürr Connect System

„The quick connection with the click!“

1. Makes assembly easier, lots of variation possible. The Dürr Connect pieces are simply plugged into each other and secured with a clip (blue, yellow or black). An o-ring on the end of the Dürr Connect parts ensures it is leak proof. The Dürr Connect parts can be used more than once – no gluing, no sealing with silicon



2. Saves space within the Dental Unit. Various different types of pieces – Elbows in 30°, 45°, 90°, T-Pieces, Y-Pieces, variable length pieces



Examples – for more see our Dürr Connect brochure P007-370-02



3. Individual connection possibilities



4. Easy and time saving for technicians when having to service the unit. Easy connection and dismantling of the units. We offer also sets where a sample of all parts is included – good for technician to keep in their car. Easy to retrofit extras such as filters, vents, rinsing unit or



5. Keeps the inside of the unit tidy



6. Easy recognition. Everything that is yellow means it is a filter and can be easily exchanged.

