


Montage-Beispiel / Mounting-Example



Montagehinweis

Bestimmungsgemäße Verwendung

DürrConnect Teile sind geeignet für dentale Anwendungen im Vakuum-, Sekret-, und Abwassersystem. Jede darüber hinausgehende Verwendung ist mit dem Hersteller abzuklären.

! Bei Verwendung der DürrConnect Teile die max. Dauergebrauchs-Temperatur und den max. Betriebsdruck beachten! Verformungen werden dadurch vermieden.

Max. Lager- u. Transport-Temperatur
- 5 °C bis + 60 °C

Max. Dauergebrauchs-Temperatur
PBT: -20 °C bis + 100 °C
PVC: -5 °C bis + 60 °C (0700-700-21/0700-800-13)

Max. Betriebsdruck
Vakuum p_{abs} : 0,5 bar / 50 kPa,
Überdruck p_g : 1,0 bar / 100 kPa

Chemikalienbeständigkeit
Beständig gegen Reinigungs- u. Desinfektionslösungen, ausgenommen Aceton

Assembly instructions

Correct usage

DürrConnect fittings are suitable for all dental plumbing purposes, including vacuum, secretion and waste water systems. Any other use must be approved by the manufacturer in advance.

! When using DürrConnect parts the max. operating temperature and the maximum operating pressure must be observed. This avoids any mal-formation taking place.

Maximum storage- and transport temperature
- 5 °C to + 60 °C

Maximum continuous operating temperature
PBT: -20 °C to + 100 °C
PVC: -5 °C to + 60 °C (0700-700-21/0700-800-13)

Maximum operating pressure
Vacuum p_{abs} : 0,5 bar / 50 kPa,
Pressure p_g : 1,0 bar / 100 kPa

Resistance to Chemicals
Resistant to most cleaning and disinfection agents with the exception of acetone.

Information

i

0700-7xx-xx
DürrConnect **klein / small**
Durchfluss / Lichte Weite
Troughflow / width inside \varnothing 20 mm

0700-8xx-xx
DürrConnect **groß / large**
Durchfluss / Lichte Weite
Troughflow / width inside \varnothing 32 mm

Übersicht DürrConnect Systeme

- 1.) DürrConnect SYSTEM Box 0700-700-50
- 2.) DürrConnect SPECIAL Box 0700-800-50
- 3.) weitere Spezial-Anschlusssteile

i Das DürrConnect Anschlusssteile Sortiment wird ständig erweitert. Informationen zu den Ersatzteilen finden Sie im Portal für autorisierte Fachhändler unter: www.duerrdental.net.

Overview DürrConnect Systems






















- 1.) DürrConnect SYSTEM Box 0700-700-50
- 2.) DürrConnect SPECIAL Box 0700-800-50
- 3.) other Special Fittings

i Our DürrConnect selection of pipe and hose connections is constantly being updated. Information on spare parts can be found on the website portal for authorised specialist dealers under: www.duerrdental.net.






















1. Anschlusssteile für / Fittings for "DürrConnect SYSTEM Box" 0700-700-50

<p>1</p> <p>5x</p> <p>0700-700-02E</p>	<p>2</p> <p>5x</p> <p>0700-700-03E</p>	<p>3</p> <p>5x</p> <p>90°</p> <p>0700-700-04E</p>	<p>4</p> <p>5x</p> <p>45°</p> <p>0700-700-05E</p>	<p>5</p> <p>5x</p> <p>30°</p> <p>0700-700-07E</p>	<p>6</p> <p>5x</p> <p>0700-700-10E</p>
<p>7</p> <p>5x</p> <p>0700-700-11E</p>	<p>8</p> <p>5x</p> <p>0700-700-13E</p>	<p>9</p> <p>2x</p> <p>0700-700-14E</p>	<p>10</p> <p>5x</p> <p>0700-700-15E</p>	<p>11</p> <p>5x</p> <p>0700-700-17E</p>	<p>12</p> <p>5x</p> <p>0700-720-10E, \varnothing10mm 0700-720-15E, \varnothing15mm 0700-720-17E, \varnothing17mm 0700-720-18E, \varnothing18mm 0700-720-20E, \varnothing20mm 0700-720-22E, \varnothing22mm 0700-720-25E, \varnothing25mm</p>
<p>13</p> <p>5x</p> <p>0700-700-20E, \varnothing32mm 0700-700-21E, \varnothing36mm</p>	<p>14</p> <p>5x</p> <p>0700-710-20E</p>	<p>15</p> <p>5x</p> <p>0700-710-25E</p>	<p>16</p> <p>9000-402-54, \varnothing20 9000-402-71, \varnothing25 9000-402-90, \varnothing30</p>	<p>17</p> <p>20x</p> <p>0700-700-01E</p>	
<p>18</p> <p>Super O-Lube Fett/Oil 2 g for O-Ring 9000-402-71</p> <p>9000474001</p>	<p>19</p> <p>Vaseline, 5 g for O-Rings 9000-402-54 + 9000-402-90</p> <p>9000-474-38</p>	<p>20</p> <p>9000-310-002</p>			

2. Anschlusssteile für / Fittings for "DürrConnect SPECIAL Box" 0700-800-50

<p>21</p>  <p>9000-402-54, ø20 9000-402-71, ø25 9000-402-90, ø30</p>	<p>22/ 23</p>  <p>20x 22: 0700-700-01E 23: 0700-800-01E</p>	<p>24</p>  <p>3x 0700-700-06E</p>	<p>25</p>  <p>4x 0700-700-16E</p>	<p>26</p>  <p>5x 0700-700-17E</p>	<p>27</p>  <p>5x 0700-700-18E</p>
<p>28</p>  <p>20x 0700-700-32E</p>	<p>29</p>  <p>5x 0700-700-33E NW16-27 0700-800-33E NW25-40</p>	<p>30</p>  <p>5x 0700-720-08E</p>	<p>31</p>  <p>5x 0700-720-30E, ø30mm</p>	<p>32</p>  <p>3x 0700-820-25E, ø25mm 0700-820-30E, ø30mm 0700-820-32E, ø32mm</p>	<p>33</p>  <p>3x 0700-800-02E</p>
<p>34</p>  <p>3x 90° 0700-800-04E</p>	<p>35</p>  <p>3x 45° 0700-800-05E</p>	<p>36</p>  <p>3x 0700-800-10E</p>	<p>37</p>  <p>3x ø 32mm ø 36mm ø 40mm 0700-800-13E</p>	<p>38</p>  <p>3x 0700-800-20E</p>	<p>39</p>  <p>5x 0700-800-21E</p>
<p>40</p>  <p>3x 0700-810-30E</p>	<p>18</p>  <p>Super O-Lube Fett/Oil 2 g for O-Ring 9000-402-71 9000474001</p>	<p>19</p>  <p>Vaseline, 5 g for O-Rings 9000-402-54 + 9000-402-90 9000-474-38</p>			

3. weitere Spezial-Anschlusssteile / other Special Fittings "DürrConnect"

<p>41</p>  <p>2x 0700-700-52E</p>	<p>42</p>  <p>10x 0700-700-54E</p>	<p>43</p>  <p>5x 0700-702-03E</p>	<p>44</p>  <p>4x 7117-420-41E</p>	<p>45</p>  <p>2x 7117-420-45E</p>	<p>50</p>  <p>3x 0700-700-08E</p>
<p>51</p>  <p>20x flach / flat 0700-700-09E</p>	<p>52</p>  <p>5x 0700-700-19E</p>	<p>53</p>  <p>0700-700-23E, ø33,5</p>	<p>54</p>  <p>5x 0700-700-26E, ø36 0700-700-27E, ø32</p>	<p>55</p>  <p>5x 0700-700-28E</p>	<p>56</p>  <p>5x 0700-700-35E</p>
<p>57</p>  <p>5x 0700-700-41E</p>	<p>58</p>  <p>2x 0700-700-43E</p>	<p>59</p>  <p>5x 0700-720-09E, ø8mm</p>	<p>60</p>  <p>5x 0700-720-23E, ø23mm</p>	<p>61</p>  <p>5x 0700-721-17E, ø17mm</p>	<p>62</p>  <p>3x 0700-810-32E</p>
<p>63</p>  <p>1x 7110-100-011</p>	<p>64</p>  <p>2x 7117-420-27E</p>	<p>65</p>  <p>3x 7128-100-03E</p>			